

## Enuresis Alarm

The alarm consists of a soft rubber mat which is placed on the bed and is connected to the bedside alarm unit. The mat is urine sensitive and triggers the battery operated alarm unit when the child wets.

For more information on our Ramsey Coote Enuresis Alarms, please visit the [Ramsey Coote Website](#).

The Alliance has three enuresis alarms for loan to our practices. These items are in great demand and there is usually a waiting list.

### Lending the Enuresis Alarm

When a request is received by the Alliance to lend the Enuresis Alarm it is referred to [Lesa Black](#). Lesa will:

1. Identify the GP who will be responsible for the equipment.
2. Produce a Lending Agreement that details the equipment to be lent.
3. Arrange for the GP to sign the Agreement.
4. Secure the \$120.00 deposit from the GP.
5. Ensure transport arrangements have been arranged by the GP for the pick-up and return of equipment.
6. Provide a copy of the Lending Agreement to the GP.

When the borrowed items are returned to the Alliance, the person receiving the goods will:

1. Check the equipment against the Lending Agreement, ensuring each item is accounted for and in the same condition as when it was lent.
2. Arrange for the return of the deposit to the GP.

### **Failure to Return Equipment**

If the borrowed equipment has not been returned by the due date the GP will first be contacted by phone, then a reminder will be sent within 3 working days.

If the equipment has not been returned within 5 working days of the due date the GP will then be invoiced for the replacement costs.

### **Terms and Conditions**

Equipment will be lent to members of the South East Alliance of General Practice (Brisbane) Ltd practicing within the Alliance's boundaries.

The Enuresis Alarm will be lent for no longer than 2 months, unless otherwise negotiated with Lesa.

A deposit for each item of equipment (totaling \$120.00) will be collected. This will be returned upon receipt of the item.

Borrowers may forfeit their deposit if any of the terms in the contract are not met.

Borrowers may be liable for the replacement cost of the equipment if the item is not returned or damaged.

The GP is entirely responsible for ensuring the Enuresis Alarm bed pad has been cleaned and sanitised.

## [Sanitising and Cleaning the Enuresis Alarm Bed Pad](#)

The manufacturer recommends the following procedure for cleaning the Bed Pad:

Practitioners immerse bed pad but not the plugs, for at least three minutes in a solution of 'Savlon' to disinfect it. Fill a basin with enough solution to cover the bed pad and lightly scrub it's surface with a nail-brush. Then lightly scrub along the electrodes with a soft toothbrush to remove lint and any other residues.

After rinsing, open-air drying is usually effective within an hour or two. Prolonged exposure to intense sunlight or heat could damage the rubber. Savlon Antiseptic Liquid, available from pharmacies, should be diluted 1:40 by volume in water. Savlon hospital concentrate, available from medical wholesalers and in larger bottles, is five times more concentrated and therefore is used in a dilution of 1:200. Other chemical agents may severely damage the pads.

It is the practitioner's responsibility to ensure the bed pad is cleaned prior to it's return to the Alliance.