

Disclaimer: Whilst every effort is made to ensure accuracy, SEA-GP (Brisbane) and acknowledged contributors do not accept any liability for any injury, loss or damage incurred by of or reliance on the information contained within this information sheet.

Background

The following information is relevant to Medical Director users only. By utilising HL7 pathology downloads and using the defined request terms, a more comprehensive patient record can be achieved.


How do I have my patient's test results appear in the diabetes record?

Diabetes Pathology test results appear in the patient's diabetes record automatically, provided that the result is sent in HL7 format and the following test names are used.

- | | | |
|----------------------------------|---------------------|--------------------|
| • INR | • Glucose | • LDL |
| • International Normalised Ratio | • SE Glucose | • LDL Cholesterol |
| • PT | • HbA1c | • HDL |
| • Prothrombin Time | • Glycosylated HB | • HDL Cholesterol |
| • Creatinine | • Glycated HB | • TG |
| • SE Creatinine | • Cholesterol | • TRIG |
| • BSL | • Total Cholesterol | • Triglycerides |
| | • SE Cholesterol | • SE Triglycerides |

Utilising the above terms and having your pathology downloaded in HL7 format will allow the graphing (as per example image) of a patient's historical results.

This can be a powerful tool to show the patient the trend in their results and facilitate patient awareness and education.



Acknowledgments: SEA-GP (Brisbane) wishes to acknowledge the valuable contributions from Dr Peter Adkins in the development of this information sheet.